

Working With Azimuth Theatre

You will have been sent a Thrive Invitation to set yourself up as part of our organization.

This app is a tool we use to better track self-assessments for Covid-19 as well as tracking Rapid Antigen Test as another layer of screening. These protocols are in place to ensure that work in person is as safe as possible.

Prior To Arrival

Complete Your Self Assessment

You can complete your self assessment before joining us in person - you can do it easily on Thrive or if you are not yet signed up by using our Google Form. More information on the Thrive app is available by scanning the first QR code.

[Thrive Info](#)




[Thrive Login](#)



[Google Form](#)



If you feel you have to answer "Yes" to any of the COVID-19 Screening questions, please contact: Azimuth Producing Team by email at producer@azimuththeatre.com or by phone at (780) 454-0583



If you have been instructed on how to complete a Rapid Antigen Test, complete it on your assigned day **prior** to joining the team in person.

If it is your first time or you are in for a one-day contract, you will be instructed on how to complete the test when you arrive.

Working With Azimuth Theatre

Upon Arrival




Mask Up

Please have a mask that covers your nose and mouth. We will always have masks available at in-person events in case you forget one!



Hand Hygiene

Wash your hands with soap and water frequently and thoroughly. You can also use hand sanitizer which will also be provided.



Temperature Checks

You will have one last temperature check before entering the space from an Azimuth Team Member.



Physical Distancing

Practice physical distancing when possible, also be aware of other folk's comfort and needs, asking and checking in on what folks' boundaries are in respectful ways.

How to test yourself with BD Veritor™ COVID-19 Rapid Antigen Screening Test

This handout explains how to use a BD Veritor™ COVID-19 rapid antigen test kit on yourself. Before you start:

If you have ANY COVID-19 symptoms or if you have been exposed to someone with COVID-19, do NOT use this test kit, and instead get tested at an assessment centre or community lab (COVID-19 testing locations)

- Make sure your kit has all the materials needed and that none of the materials are expired or damaged.
- Store the kit and its content at room temperature (not in the fridge or freezer) and away from direct sunlight.
- Find a clean area free of foods, drinks, or clutter, and easy to clean in the event of a spill.
- During the test, do not eat, drink, smoke, vape, put in contact lenses, put on make-up, or touch your face in any way.
- This test is for your own use only, unless you are a parent/guardian whose child is being tested.

1 Blow your nose and then wash your hands for at least 20 seconds.


2 Assemble the materials in your kit. They include:




One test strip



One sterile nasal swab



Tube with reagent fluid*
(*See Note on page 2)




3 Swab both nostrils

- Remove the swab from its package, keeping the swab tip away from your hands or any other surfaces.
- Tilt your head back and insert the swab straight back (not up) into one nostril. Stop when you meet resistance or when you reach 2.5 cm (distance is shorter in children < 12 years old).
- Turn the swab against the inside of your nose 5-10 times.
- Remove the swab and repeat same steps in the other nostril.


4 Prepare the Sample



• While making sure that the swab tip never touches any surface or your hands, remove cap from the tube and insert the nasal swab



• Swirl the swab in the liquid for 15 seconds, pressing it against the sides of the tube.



• Use one hand to pinch the tube and pull up the swab with the other hand, bringing the swab through the pinched tube.




• Throw the swab in the garbage.
• Firmly close the tube with the dispensing cap (attached to the tube).


How to test yourself with BD Veritor™ COVID-19 Rapid Antigen Screening Test

5 Perform the test.

- Remove test strip from its package and place it face up on a clean and flat surface.
- Turn the tube upside down over the square sample well on the test strip and gently squeeze the tube to put 3 drops into the sample well (DO NOT put drops in the result window).



Scan Code for video support




Result window

6 Wait 15 minutes before reading your results.


- Your result may not be correct if read in less than 15 minutes or more than 20 minutes after putting the drops into the sample well.
- Using a marker, write your initials and date on the test.



7 Read your results. It does not matter which line appeared first or how faint the lines look.



• If only the control line (C) is present, the result is negative.



• If both control line (C) and test line (T) are present, the result is positive.



or



• If the control line (C) is not present, result is not valid (i.e., the test is not working).

9 Dispose of test materials in the garbage and wash your hands.

10 Take a picture of the test result and upload to thrive

8 Act on your result quickly.

• If your result is negative, it is less likely that you have COVID-19. You may resume your daily activities while continuing to follow all public health advice.

• If your result is positive, you may have COVID-19 and should get a laboratory test at an assessment centre or community lab as soon as possible to confirm your result. Immediately self-isolate until you receive your laboratory test result.

• If your test is not valid, repeat the test using a new swab and test device. You will need to swab your nose again according to the instructions.

*Note: The liquid (reagent) in the tubes can irritate (bother) your eyes and skin. The liquid is toxic (dangerous) if swallowed. If the liquid accidentally touches your skin, wash it with plenty of water. If it touches your eyes or mouth, rinse with plenty of water for 15 minutes and contact your health care provider.